


Product datasheet for **GA118338**

Human WDR38 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR38 CRISPRa kit - CRISPR gene activation of human WD repeat domain 38
Symbol:	WDR38
Components:	GA118338G1, WDR38 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118338G2, WDR38 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118338G3, WDR38 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5JTN6
Synonyms:	OTTHUMP00000022116; WD repeat domain 38
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)