

Product datasheet for **GA118355**

Human SRRD activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SRRD CRISPRa kit - CRISPR gene activation of human SRR1 domain containing
Symbol:	SRRD
Components:	GA118355G1, SRRD gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118355G2, SRRD gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118355G3, SRRD gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UH36</u>
Synonyms:	HC/HCC; SRR1L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)