

## Product datasheet for **GA118367**

### Human IFITM10 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | IFITM10 CRISPRa kit - CRISPR gene activation of human interferon induced transmembrane protein 10                                                                                                                                                                                                                |
| Symbol:       | IFITM10                                                                                                                                                                                                                                                                                                          |
| Components:   | GA118367G1, IFITM10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA118367G2, IFITM10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA118367G3, IFITM10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>A6NMD0</u>                                                                                                                                                                                                                                                                                                    |
| Synonyms:     | DSPA3                                                                                                                                                                                                                                                                                                            |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)