

Product datasheet for **GA118389**

Human SCGB1D4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SCGB1D4 CRISPRa kit - CRISPR gene activation of human secretoglobin family 1D member 4
Symbol:	SCGB1D4
Components:	GA118389G1, SCGB1D4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118389G2, SCGB1D4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118389G3, SCGB1D4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6XE38</u>
Synonyms:	IIS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)