

Product datasheet for **GA118396**

Human TAF_{A1} activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF _{A1} CRISPRa kit - CRISPR gene activation of human TAF _A chemokine like family member 1
Symbol:	TAF _{A1}
Components:	GA118396G1, TAF _{A1} gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118396G2, TAF _{A1} gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118396G3, TAF _{A1} gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z5A9</u>
Synonyms:	TAF _A -1; TAF _{A1}
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)