

Product datasheet for **GA118425**

Human MXRA7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MXRA7 CRISPRa kit - CRISPR gene activation of human matrix remodeling associated 7
Symbol:	MXRA7
Components:	GA118425G1, MXRA7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118425G2, MXRA7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118425G3, MXRA7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P84157</u>
Synonyms:	PS1TP1; TMAP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)