


Product datasheet for **GA118426**

Human IFIT1B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IFIT1B CRISPRa kit - CRISPR gene activation of human interferon induced protein with tetratricopeptide repeats 1B
Symbol:	IFIT1B
Components:	GA118426G1, IFIT1B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118426G2, IFIT1B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118426G3, IFIT1B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T764</u>
Synonyms:	bA149I23.6; IFIT1L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)