

Product datasheet for GA118429

Human TMEM41B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM41B CRISPRa kit - CRISPR gene activation of human transmembrane protein 41B
Symbol:	TMEM41B
Components:	GA118429G1, TMEM41B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118429G2, TMEM41B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118429G3, TMEM41B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5BJD5
Synonyms:	KIAA0033; MGC33897; MGC57262
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)