

## Product datasheet for **GA118436**

### Human Histone H3.3C (H3F3C) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	H3F3C CRISPRa kit - CRISPR gene activation of human H3 histone family member 3C
Symbol:	H3-5
UniProt ID:	<u><a href="#">Q6NXT2</a></u>
Synonyms:	H3.5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)