

Product datasheet for **GA118447**

Human GCN2 (EIF2AK4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF2AK4 CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 2 alpha kinase 4
Symbol:	EIF2AK4
UniProt ID:	<u>Q9P2K8</u>
Synonyms:	GCN2; PVOD2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)