

Product datasheet for **GA118491**

Human TMEM151B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM151B CRISPRa kit - CRISPR gene activation of human transmembrane protein 151B
Symbol:	TMEM151B
Components:	GA118491G1, TMEM151B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118491G2, TMEM151B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118491G3, TMEM151B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8IW70
Synonyms:	bA444E17.5; C6orf137; TMEM193
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)