

Product datasheet for **GA118498**

Human SPDYE3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPDYE3 CRISPRa kit - CRISPR gene activation of human speedy/RINGO cell cycle regulator family member E3
Symbol:	SPDYE3
Components:	GA118498G1, SPDYE3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118498G2, SPDYE3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118498G3, SPDYE3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A6NKU9</u>
Synonyms:	MGC57359; MGC119294
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)