


Product datasheet for **GA118516**

Human CT45A3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CT45A3 CRISPRa kit - CRISPR gene activation of human cancer/testis antigen family 45 member A3
Symbol:	CT45A3
Components:	GA118516G1, CT45A3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118516G2, CT45A3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118516G3, CT45A3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NHU0</u>
Synonyms:	CT45-3; CT45-4; CT45.3; CT45.4; CT45A4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)