

Product datasheet for **GA118519**

Human **PGAM4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	PGAM4 CRISPRa kit - CRISPR gene activation of human phosphoglycerate mutase family member 4
Symbol:	PGAM4
UniProt ID:	<u>Q8N0Y7</u>
Synonyms:	dj1000K24.1; PGAM-B; PGAM1; PGAM3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)