


Product datasheet for **GA118550**

Human LMOD2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LMOD2 CRISPRa kit - CRISPR gene activation of human leiomodin 2
Symbol:	LMOD2
Components:	GA118550G1, LMOD2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118550G2, LMOD2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118550G3, LMOD2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P5Q4</u>
Synonyms:	C-LMOD; CLMOD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)