

Product datasheet for **GA118555**

Human ARHGEF35 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARHGEF35 CRISPRa kit - CRISPR gene activation of human Rho guanine nucleotide exchange factor 35
Symbol:	ARHGEF35
Components:	GA118555G1, ARHGEF35 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118555G2, ARHGEF35 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118555G3, ARHGEF35 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A5YM69
Synonyms:	ARHGEF5L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)