

Product datasheet for **GA118564**

Human KPRP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KPRP CRISPRa kit - CRISPR gene activation of human keratinocyte proline rich protein
Symbol:	KPRP
Components:	GA118564G1, KPRP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118564G2, KPRP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118564G3, KPRP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T749</u>
Synonyms:	C1orf45
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)