


## Product datasheet for **GA118566**

### Human **LRRC37A2** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	LRRC37A2 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 37 member A2
Symbol:	LRRC37A2
Components:	GA118566G1, LRRC37A2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118566G2, LRRC37A2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118566G3, LRRC37A2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">A6NM11</a></u>
Synonyms:	LRRC37
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)