

Product datasheet for **GA118569**

Human **LRRC18** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC18 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 18
Symbol:	LRRC18
Components:	GA118569G1, LRRC18 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118569G2, LRRC18 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118569G3, LRRC18 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N456</u>
Synonyms:	UNQ933; UNQ9338; VKGE9338
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)