

Product datasheet for **GA118573**

Human F8A3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	F8A3 CRISPRa kit - CRISPR gene activation of human coagulation factor VIII associated 3
Symbol:	F8A3
Components:	GA118573G1, F8A3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118573G2, F8A3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118573G3, F8A3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P23610</u>
Synonyms:	coagulation factor VIII-associated (intronic transcript) 3; OTTHUMP00000032188
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)