

Product datasheet for **GA118579**

Human **CCDC73** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CCDC73 CRISPRa kit - CRISPR gene activation of human coiled-coil domain containing 73
Symbol:	CCDC73
Components:	GA118579G1, CCDC73 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118579G2, CCDC73 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118579G3, CCDC73 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZRK6</u>
Synonyms:	NY-SAR-79
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)