

Product datasheet for **GA118580**

Human EID3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EID3 CRISPRa kit - CRISPR gene activation of human EP300 interacting inhibitor of differentiation 3
Symbol:	EID3
UniProt ID:	Q8N140
Synonyms:	NS4EB; NSE4B; NSMCE4B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)