

Product datasheet for **GA118587**

Human CHAC2 activation kit by CRISPRa

Product data:

Product Type: CRISPRa Kits

Description: CHAC2 CRISPRa kit - CRISPR gene activation of human ChaC cation transport regulator homolog 2

Symbol: CHAC2


Components: GA118587G1, CHAC2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)
GA118587G2, CHAC2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)
GA118587G3, CHAC2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)
1 CRISPRa-Enhancer vector, SKU GE100056
1 CRISPRa scramble vector, SKU GE100077

UniProt ID: Q8WUX2

Synonyms: ChaC, cation transport regulator-like 2; ChaC, cation transport regulator homolog 2 (E. coli)

Format: 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)