

Product datasheet for **GA118597**

Human MS4A13 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MS4A13 CRISPRa kit - CRISPR gene activation of human membrane spanning 4-domains A13
Symbol:	MS4A13
Components:	GA118597G1, MS4A13 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118597G2, MS4A13 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118597G3, MS4A13 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5J8X5
Synonyms:	membrane-spanning 4-domains, subfamily A, member 13; testis-expressed transmembrane-4 protein
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)