

Product datasheet for **GA118639**

Human CYB5RL activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CYB5RL CRISPRa kit - CRISPR gene activation of human cytochrome b5 reductase like
Symbol:	CYB5RL
Components:	GA118639G1, CYB5RL gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118639G2, CYB5RL gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118639G3, CYB5RL gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6IPT4</u>
Synonyms:	FLJ16295
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)