

Product datasheet for **GA118646**

Human HIGD1C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIGD1C CRISPRa kit - CRISPR gene activation of human HIG1 hypoxia inducible domain family member 1C
Symbol:	HIGD1C
Components:	GA118646G1, HIGD1C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118646G2, HIGD1C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118646G3, HIGD1C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A8MV81</u>
Synonyms:	GM921
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)