


Product datasheet for **GA118647**

Human SERINC4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SERINC4 CRISPRa kit - CRISPR gene activation of human serine incorporator 4
Symbol:	SERINC4
Components:	GA118647G1, SERINC4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118647G2, SERINC4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118647G3, SERINC4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A6NH21</u>
Synonyms:	FLJ40363
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)