


Product datasheet for **GA118649**

Human ZNF704 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF704 CRISPRa kit - CRISPR gene activation of human zinc finger protein 704
Symbol:	ZNF704
Components:	GA118649G1, ZNF704 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118649G2, ZNF704 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118649G3, ZNF704 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZNC4</u>
Synonyms:	Gig1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)