

Product datasheet for **GA118678**

Human TMEM232 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM232 CRISPRa kit - CRISPR gene activation of human transmembrane protein 232
Symbol:	TMEM232
Components:	GA118678G1, TMEM232 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118678G2, TMEM232 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118678G3, TMEM232 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	C9JQI7
Synonyms:	DKFZp686G0585; FLJ35903; FLJ43080
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)