

Product datasheet for **GA118693**

Human SPINK9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPINK9 CRISPRa kit - CRISPR gene activation of human serine peptidase inhibitor, Kazal type 9
Symbol:	SPINK9
Components:	GA118693G1, SPINK9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118693G2, SPINK9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118693G3, SPINK9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5DT21</u>
Synonyms:	LEKT12
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)