

## Product datasheet for **GA118695**

### Human LIPN activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	LIPN CRISPRa kit - CRISPR gene activation of human lipase family member N
Symbol:	LIPN
Components:	GA118695G1, LIPN gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118695G2, LIPN gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118695G3, LIPN gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5VXI9</u>
Synonyms:	ARCI8; bA186O14.3; LI4; LIPL4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)