

Product datasheet for **GA118708**

Human **PGA3** activation kit by **CRISPRa**

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PGA3 CRISPRa kit - CRISPR gene activation of human pepsinogen 3, group I (pepsinogen A) |
| Symbol: | PGA3 |
| UniProt ID: | <u>P00790</u> |
| Synonyms: | PGA4; PGA5 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)