

Product datasheet for **GA118713**

Human CBLN3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CBLN3 CRISPRa kit - CRISPR gene activation of human cerebellin 3 precursor
Symbol:	CBLN3
Components:	GA118713G1, CBLN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118713G2, CBLN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118713G3, CBLN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6UW01</u>
Synonyms:	PRO1486
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)