

Product datasheet for **GA118747**

Human **GAGE13** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	GAGE13 CRISPRa kit - CRISPR gene activation of human G antigen 13
Symbol:	GAGE13
Components:	GA118747G1, GAGE13 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118747G2, GAGE13 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118747G3, GAGE13 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q4V321</u>
Synonyms:	G antigen 12D; MGC163246; OTTHUMP00000025843
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)