

Product datasheet for **GA118754**

Human LINGO3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LINGO3 CRISPRa kit - CRISPR gene activation of human leucine rich repeat and Ig domain containing 3
Symbol:	LINGO3
Components:	GA118754G1, LINGO3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118754G2, LINGO3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118754G3, LINGO3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P0C6S8</u>
Synonyms:	LERN2; LRRN6B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)