

Product datasheet for **GA118757**

Human TTMA (TMEM200C) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM200C CRISPRa kit – CRISPR gene activation of human transmembrane protein 200C
Symbol:	TMEM200C
Components:	GA118757G1, TTMA gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118757G2, TTMA gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118757G3, TTMA gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A6NKL6
Synonyms:	TTMA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

