

Product datasheet for **GA118775**

Human SPINK8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPINK8 CRISPRa kit - CRISPR gene activation of human serine peptidase inhibitor, Kazal type 8 (putative)
Symbol:	SPINK8
Components:	GA118775G1, SPINK8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118775G2, SPINK8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118775G3, SPINK8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P0C7L1</u>
Synonyms:	C630041L24Rik; serine peptidase inhibitor, Kazal type 8; serine protease inhibitor Kazal type 8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)