

Product datasheet for **GA118783**

Human SBK2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SBK2 CRISPRa kit - CRISPR gene activation of human SH3 domain binding kinase family member 2
Symbol:	SBK2
Components:	GA118783G1, SBK2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118783G2, SBK2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118783G3, SBK2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	SGK069
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)