


Product datasheet for **GA118784**

Human TMEM90A (SYNDIG1L) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SYNDIG1L CRISPRa kit - CRISPR gene activation of human synapse differentiation inducing 1 like
Symbol:	SYNDIG1L
Components:	GA118784G1, TMEM90A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118784G2, TMEM90A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118784G3, TMEM90A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A6NDD5
Synonyms:	CAPUCIN; DSPC1; IFITMD4; TMEM90A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)