

Product datasheet for **GA118809**

Human PPP1R3G activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PPP1R3G CRISPRa kit - CRISPR gene activation of human protein phosphatase 1 regulatory subunit 3G |
| Symbol: | PPP1R3G |
| UniProt ID: | <u>B7ZBB8</u> |
| Synonyms: | 1600032L17Rik; protein phosphatase 1, regulatory (inhibitor) subunit 3G |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:

