

Product datasheet for **GA118836**

Human KIAA0895L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KIAA0895L CRISPRa kit - CRISPR gene activation of human KIAA0895 like
Symbol:	KIAA0895L
Components:	GA118836G1, KIAA0895L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118836G2, KIAA0895L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118836G3, KIAA0895L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q68EN5</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

