

Product datasheet for **GA118847**

Human PPIAL4A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPIAL4A CRISPRa kit - CRISPR gene activation of human peptidylprolyl isomerase A like 4A
Symbol:	PPIAL4A
UniProt ID:	Q9Y536
Synonyms:	COAS-2; COAS2; PPIAL4; PPIAL4A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)