

## Product datasheet for **GA118850**

### Human TP53TG3C activation kit by CRISPRa

#### Product data:

<b>Product Type:</b>	CRISPRa Kits
<b>Description:</b>	TP53TG3C CRISPRa kit - CRISPR gene activation of human TP53 target 3C
<b>Symbol:</b>	TP53TG3C
<b>Components:</b>	GA118850G1, TP53TG3C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118850G2, TP53TG3C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118850G3, TP53TG3C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
<b>UniProt ID:</b>	<u><a href="#">Q9ULZ0</a></u>
<b>Synonyms:</b>	TP53TG3; TP53TG3E; TP53TG3F
<b>Format:</b>	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

