

## Product datasheet for **GA118851**

### Human **FAM23A (TMEM236)** activation kit by **CRISPRa**

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | TMEM236 CRISPRa kit – CRISPR gene activation of human transmembrane protein 236                                                                                                                                                                                                                               |
| <b>Symbol:</b>       | TMEM236                                                                                                                                                                                                                                                                                                       |
| <b>Components:</b>   | GA118851G1, FAM23A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA118851G2, FAM23A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA118851G3, FAM23A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <a href="#">Q5W0B7</a>                                                                                                                                                                                                                                                                                        |
| <b>Synonyms:</b>     | bA16O1.2; bA162I21.2; FAM23A; FAM23B                                                                                                                                                                                                                                                                          |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

