

## Product datasheet for **GA118853**

### Human PPIAL4C activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PPIAL4C CRISPRa kit - CRISPR gene activation of human peptidylprolyl isomerase A like 4C
Symbol:	PPIAL4C
UniProt ID:	<a href="https://www.uniprot.org/entry/A0A0B4J2A2">A0A0B4J2A2</a>
Synonyms:	COAS-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)