

Product datasheet for **GA118861**

Human TMEM183BP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM183B CRISPRa kit - CRISPR gene activation of human transmembrane protein 183B
Symbol:	TMEM183BP
Components:	GA118861G1, TMEM183BP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118861G2, TMEM183BP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118861G3, TMEM183BP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q1AE95
Synonyms:	C1orf37-DUP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)