

Product datasheet for **GA118877**

Human PCP4L1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCP4L1 CRISPRa kit - CRISPR gene activation of human Purkinje cell protein 4 like 1
Symbol:	PCP4L1
Components:	GA118877G1, PCP4L1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118877G2, PCP4L1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118877G3, PCP4L1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A6NKN8</u>
Synonyms:	IQM1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)