

Product datasheet for **GA118913**

Human TSPY3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSPY3 CRISPRa kit - CRISPR gene activation of human testis specific protein Y-linked 3
Symbol:	TSPY3
Components:	GA118913G1, TSPY3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118913G2, TSPY3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118913G3, TSPY3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q01534</u>
Synonyms:	CT78
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)