

Product datasheet for **GA118918**

Human TMEM191B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM191B CRISPRa kit – CRISPR gene activation of human transmembrane protein 191B
Symbol:	TMEM191B
Components:	GA118918G1, TMEM191B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118918G2, TMEM191B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118918G3, TMEM191B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P0C7N4</u>
Synonyms:	AC023490.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

