

Product datasheet for **GA118941**

Human TSPY8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSPY8 CRISPRa kit – CRISPR gene activation of human testis specific protein Y-linked 8
Symbol:	TSPY8
Components:	GA118941G1, TSPY8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA118941G2, TSPY8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA118941G3, TSPY8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P0CW00</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

